



nSTRIDE[®] Autologous Protein Solution (APS) Kit

Once OA Pain Starts, It's Hard to Stop.



ZIMMER BIOMET
Your progress. Our promise.[®]



nSTRIDE® Autologous Protein Solution (APS) Kit

The Knee OA Treatment Gap

It is well recognized that a treatment gap exists when considering solutions for osteoarthritic pain.¹² Between conservative care and the more invasive implant solutions there are few well researched and clinically proven products that can relieve the patient's pain and potentially delay the need for an implant.

Autologous anti-inflammatory Injection for the Treatment of Knee Osteoarthritis

Once OA pain starts it is hard to stop. The nSTRIDE APS Kit is designed to produce a novel therapy to treat pain and slow the progression of cartilage degradation and destruction in the knee. The nSTRIDE APS Kit is a point of care cell-concentration system which concentrates anti-inflammatory cytokines and anabolic growth factors to significantly decrease pain and promote cartilage health.

Extensive Research Program and clinical results

nSTRIDE Autologous Protein Solution (APS) was developed after a multi-year research program which focused on understanding the osteoarthritic disease process in the knee and understanding the Mode of Action that an autologous blood based product could have. With the nSTRIDE Kit an Autologous Anti-Inflammatory (AAI) solution is produced from the patient's own blood and this is injected intra articularly into the knee.

nSTRIDE APS has been shown to

- Significantly Reduces Pain Associated with Knee OA up to 2 years^{1-3,13}
- Significantly Improves function in the Knee Joint associated with OA¹⁻³
- To be effective for patients with Kellgren and Lawrence stage 2 and 3 following a single injection^{1-3,13}
- 70% Improvement in Knee Pain at 2 years following a Single Injection¹³

“There is great need for a safe, effective, and cost effective treatment option for patients with moderate to severe osteoarthritis that enjoys high patient acceptance”.¹²




Figure 1:
 In an osteoarthritic knee, an increase in inflammatory cytokines results in cartilage degeneration and knee pain. The inflammatory proteins IL-1 and TNF α attack the cartilage. These inflammatory proteins must be stopped simultaneously to decrease pain and slow cartilage degeneration.⁵

Understanding the Mode of Action of Knee Osteoarthritis

In an osteoarthritic knee, an increase in inflammatory cytokines results in cartilage degeneration and knee pain. The inflammatory proteins IL-1 and TNF α attack the cartilage. These inflammatory proteins must be stopped simultaneously to decrease pain and slow cartilage degeneration⁵.

Understanding the Mode of Action

In pre-clinical studies, nSTRIDE APS has been shown to:

- Inhibit catabolic enzyme production from chondrocytes stimulated with IL-1 β and TNF α ^{8^}
- Inhibit inflammatory cytokine production from IL-1 β -stimulated macrophages^{9^}
- Inhibit catabolic destruction of cartilage tissue^{4^*}
- Protect cartilage in a meniscal-tear model^{10^*}
- Reduce pain in large animals with naturally occurring OA^{11^*}
- Stimulates Cartilage Cell Proliferation^{4^*}




Figure 2: Cytokine Imbalance in OA¹⁴








Figure 3: The nSTRIDE APS output introduces high levels of IL-1ra, sIL-1R, sTNF-RI, and sTNF-RII that block the inflammatory cytokines IL-1 and TNF α^7 which may slow cartilage degeneration.^{4**} While balance is being restored to the knee, anabolic growth factors (IGF-1 and TGF- β 1)⁷ are also introduced for beneficial cartilage health.⁶

How does nSTRIDE APS work?

Autologous Anti-Inflammatory Therapy

The nSTRIDE APS Kit processes a patient's own blood at the point-of-care to concentrate white blood cells, platelets, and plasma proteins into a small volume of plasma. The output contains the anti-inflammatory cytokines IL-1ra, sIL-1R, sTNF-RI, and sTNF-RII in concentrations well above what is found in native whole blood.⁶ In addition to the anti-inflammatory cytokines, anabolic cytokines for cartilage including IGF-1 and TGF- α are also concentrated to levels well in excess of that found in whole blood⁷.

The APS produced by the nSTRIDE APS Kit reduces prevalent inflammatory cytokine activity which is upregulated in osteoarthritis. These cytokines can cause both pain and cartilage degeneration.⁶ Inflammation is implicated in both the pain and the cartilage matrix breakdown in osteoarthritic joints. The cytokines IL-1 and TNF α are the key pro-inflammatory and catabolic targets; however, they need to be inhibited simultaneously.⁵ There are several naturally occurring inhibitors of IL-1 and TNF α , including IL-1ra, sIL-1R, sTNF-RI, and sTNF-RII. Additionally, while decreased activity of IL-1 and TNF α will decrease inflammation and slow the progression of cartilage degradation, anabolic cytokines could also assist to stimulate cartilage matrix synthesis.^{4**}

Therefore, the ideal therapy to target pain and cartilage matrix degradation caused by inflammation will include several different anti-inflammatory cytokines and anabolic cytokines together. Autologous Protein Solution, produced using the nSTRIDE Autologous Protein Solution Kit, contains all the anti-inflammatory proteins mentioned above, with an added benefit of anabolic cytokines as well.⁷




Figure 4: nSTRIDE APS Addresses Cytokine Imbalance in OA¹⁴

Long lasting results




Figure 5: WOMAC percent improvement with nSTRIDE APS. 69.7% improvement, ($p < 0.0001$) 24 months post-injection[#]

70% Improvement in Knee Pain
at 2 years following a Single
Injection^{13 #}

Clinical Results

In a series of clinical studies the autologous anti-inflammatory solution produced by the nSTRIDE APS kit has been shown to have positive effects on knee pain and function in patients with Kellgren-Lawrence stage 2 and 3. Further key findings were:

- Significantly Reduces Pain Associated with Knee OA up to 2 years following a single injection^{1-3,13} (Figure 1)
- Significant pain reduction using leukocyte-containing APS compared to baseline^{1,2}
- IL-1ra concentration correlated to white blood cell content in APS¹
- L-1β concentration did not increase with white blood cell content¹
- The ratio of IL-1ra to IL-1β in APS was significantly correlated with improved WOMAC pain scores at six months post-injection¹
- 72.7% of subjects were OMERACT-OARSI high responders six months post-injection²
- Significant percent improvement in pain over saline injection in double blinded pilot study^{3,13}



References

1. King W, van der Weegen W, Van Drumpt R, Soons H, Toler K, Woodell-May J, "White bloodcell concentration correlates with increased concentration of IL-1ra and improvement in WOMAC pain scores in an open-label safety study of autologous protein solution." Journal of Experimental Orthopaedics. 2016;3:9.
 2. van Drumpt RA, van der Weegen W, King WJ, Toler K, Macenski M. Safety and treatment effectiveness of a single autologous protein solution injection in patients with knee osteoarthritis. Accepted for publication 2016.
 3. Kon E, Engebretsen L, Peter Verdonk P, Nehrer S and Filardo G. "Clinical Outcomes of Knee Osteoarthritis Treated with an Autologous Protein Solution. A 1-year Pilot Double-Blinded Randomized Control Trial. American Journal of Sports Medicine, Oct. 2017.
 4. Matuska A, O'Shaughnessey KM, King WJ, Woodell-May JE. Autologous solution protects bovine cartilage explants from IL-1 β and TNF α induced cartilage degradation. J Orthop Res ,2013;31(12):1929-35.
 5. Abbas AK, Lichtman A. Cytokines. Cellular and Molecular Immunology. 5 ed. Philadelphia: Saunders; 2003. p. 243-74.
 6. Goldring SR, Goldring MB. The role of cytokines in cartilage matrix degeneration in osteoarthritis. Clin Orthop Relat Res 2004 Oct;(427 Suppl):S27-S36.
 7. O'Shaughnessey K, Matuska A, Hoepfner J, et al. Autologous protein solution prepared from the blood of osteoarthritic patients contains an enhanced profile of anti-inflammatory cytokines and anabolic growth factors. J Orthop Res 2014 Oct;32(10):1349-55.
 8. Woodell-May J, Matuska A, Oyster M, et al. Autologous protein solution inhibits MMP-13 production by IL-1beta and TNFalpha-stimulated human articular chondrocytes. J Orthop Res 2011 Sep 15;29:1320-6.
 9. O'Shaughnessey K, Panitch A, Woodell-May JE. Blood-derived anti-inflammatory protein solution blocks the effect of IL-1 β on human macrophages in vitro. Inflamm Res 2011 Oct;60(10):929-36.
 10. King W., Bendele A., Marohl T., Woodell-May J., "Human blood-based anti-inflammatory solution inhibits osteoarthritis progression in a meniscal-tear rat study.," Journal of Orthopaedic Research, 35(10):2260-2268., 2017.
 11. Bertone AL, Ishihara A, Zekas LJ, et al. Evaluation of a single intra-articular injection of autologous protein solution for treatment of osteoarthritis in horses. Am J Vet Res 2014 Feb 1;75(2):141-51.
 12. Nicholas J. London, Larry E. Miller, Jon E. Block Med Hypotheses. 2011 76(2011) 887-892
 13. Kon E, Engebretsen L, Peter Verdonk P, Nehrer S and Filardo G. "Two-year Clinical Outcomes of An Autologous Protein Solution Injection For Knee Osteoarthritis." ICRS 14th World Congress, presented, 2018.
 14. Goldring MB. The role of the chondrocyte in osteoarthritis. Arthritis Rheum. 2000;43(9):1916-1926
 15. Hix J, Klaassen M, Foreman R, Cullen E, Toler K, King W, Woodell-May J, " An Autologous Anti-Inflammatory Protein Solution Yielded a Favorable Safety Profile and Significant Pain Relief in an Open-Label Pilot Study of Patients with Osteoarthritis," BioResearch Open Access,6(1):151-158. 2017.
- ^ Cell culture assays are not necessarily indicative of clinical outcomes.
* Animal studies are not necessarily indicative of clinical outcomes
As measured by WOMAC pain scores reported by patients continuing follow-up through 2 years (n = 22)

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. This material is intended for health care professionals and the Zimmer Biomet sales force. Distribution to any other recipient is prohibited. For complete product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert and www.zimmerbiomet.com. Not intended for surgeons practicing medicine in France. Zimmer Biomet does not practice medicine. The treating surgeon is responsible for determining the appropriate treatment, techniques(s), and product(s) for each individual patient.

©2018 Zimmer Biomet. All rights reserved.



2062.1-EMEA-en-REV0518

Responsible Manufacturer

Zimmer Biomet
P.O. Box 587
56 E. Bell Drive
Warsaw, Indiana 46581-0587
USA
www.zimmerbiomet.com

Authorized Representative

Zimmer Biomet UK
Waterton Industrial Estate
Bridgend, South Wales
CF31 3XA
UK

